

# MEASUREMENT OF PRESSURE CHANGE WITH DENSITY CHANGE

MED 16.03



## Material:

Item-no.	Qty.	Description
C3020-6D	1	Erlenmeyer flask glass, 250 ml, SB 29
C7320-4B	1	Stopper silicone, 26/32/30 mm, 1 hole, for SB 29
C6230-1B	1	Stopcock 3 way, PP
C7445-7G	1	Hose, rubber, D=7/10 mm, L=100 cm
C7445-5S	1	Hose, silicone, D=5/8 mm, L=100 cm
C6095-1K	1	Tubing connector for di= 7-10 mm
C6100-2A	1	Syringe, 120ml, plastics
DE722-2D	1	Manometer, differential, "inno"
P3120-4A	1	L-shaped assembly platform

# MEASUREMENT OF PRESSURE CHANGE WITH DENSITY CHANGE

MED 16.03

## Purpose

Measuring the increase of pressure increase with increasing density.

## Preparation

- the three-way valve is inserted into the stopper as shown on the image
- insert the stopper into the Erlenmeyer flask
- the outlet on the side of the valve is connected to the Manometer via a silicone hose
- the piston of the syringe is pulled outwards and connect it with the upper end of the valve via another silicone hose
- the three-way valve is set so that there is a connection between the piston and the measuring device (measuring position)
- turn on the Manometer and set a measuring range of „1000 hPa“, afterwards tare (set to 0)



## Experiment

The three-way valve is set so that there is a connection between the Erlenmeyer flask and the piston. Push air with the piston into the Erlenmeyer flask.

Afterwards bring the three-way valve back into „measuring position“ and read off the result on the Manometer.

## Result

After turning to the measuring position while maintaining the piston pressure, the increased pressure is determined - caused by the larger number of impacting particles per unit area.

## Note

If you also want to establish a relationship between the increase in density and the increase in pressure, special attention must be paid to a constant temperature.