

MEASURING VOLUME OF LIQUIDS

MED 01.03



Material:

Item-no.	Qty.	Description
C1580-1B	1	Measuring cylinder glass, 100 ml
DG110-1B	1	Measuring beaker, plastic, 1000 ml
	1	Coloured Water

MEASURING VOLUME OF LIQUIDS

MED 01.03

Purpose

Determination of liquid volumes - practice in scale reading

Experiment

Pour 150 ml of water into the measuring beaker.

Afterwards pour 20 ml, then 78 ml and finally 100 ml water into the measuring cylinder.

The correct reading on the graduation of the measuring cylinder must be ensured.